

UV-VIS SPECTROSCOPY



Spectrophotometer

INTRODUCTION

- ***Spectroscopy*** is the branch of science which deals with interaction of electromagnetic radiation with materials. In other words it is an analytical method for qualitative and quantitative analysis by use of light.
- ***The Lambert Beer law*** in 1852 made the basis for the quantitative evaluation of absorption measurements.²



1.SPECTROSCOPY

- When an Electromagnetic radiation is incident on a matter, phenomena like reflection, transmission, absorption ,are occurring.⁴
- Spectroscopy is the study of interaction of electromagnetic radiation with matter based on the Bohr-Einstein frequency relationship $E=h\nu$, here h is the proportionality constant called Planck's constant ($6.626 \times 10^{-34} \text{ J s}$) and ν is frequency.



- Measurement of radiation intensity as a function of wavelength is described by ***spectroscopy***, as shown in figure 2.




FIG.2.Spectroscopy graph.

- All forms of ***spectroscopy*** use part of the electromagnetic radiation to give us information about the materials.

1.1 SPECTRUM

- The *spectrum* is formed by electromagnetic waves and the wavelength is varies. See figure 3.




FIG.3.Electromagnetic spectrum.

- When a narrow beam of light is allowed to pass through a prism/grating, it is dispersed into seven colors from red to violet and the band is called *Spectrum*. See figure 4.




FIG.4.Glase prism dispersion

1.2 PRINCIPLE

Basic principle of *spectroscopy* is the *Beer-Lambert's law*.²

1.2.1 BEER LAW


- *Beer's law* stated that **absorbance** is proportional to the **concentrations** of the material sample.

1.2.2 LAMBERT LAW

- *Lambert's law* stated that **absorbance** of a material is directly proportional to its **thickness** (path length).



- The modern derivation of the Beer–Lambert law combines the two laws and correlates the absorbance to both the concentrations and the thickness of the material.



$$A = \epsilon b c$$

A = Absorbance = $-\log_{10} I / I_0$

b = Optical path length (cm)

c = Solution Concentration (M)

ϵ = Molar Absorptivity $1/(\text{cm} \cdot \text{M})$

(measure of how well a compound absorbs a given wavelength of light)

I_0 = Incident beam
intensity

I = Transmitted beam
intensity

FIG 5. Beer–Lambert law.

2. UV-VISIBLE SPECTROSCOPY

- ***Ultraviolet-visible spectrum*** can be generated when ultraviolet light and visible light(200-900nm) are absorbed by materials. The spectrum can be used to analyze the composition and the structure of the materials. For a particular wavelength in the ultraviolet-visible ranges, the absorption degree is proportional to the components of the materials. Therefore, the characteristics of the materials are quantitatively reflected by the spectrum, which changes with the wave-length.¹
- ***Ultraviolet-visible spectrum*** consists of an absorption spectrum. An absorption spectrum gives information about the *molar absorptivity*, *concentration of the sample*, *optical bath length*. See figure6, in previous slid.



2.1 INSTRUMENTATION

- 2.1.1 SOURCE of LIGHT.
- 2.1.2 MONOCHROMATOR.
- 2.1.3 SMPLE SOLIOTION in CUVETTE.
- 2.1.4 PHOTO DETECTOR.
- 2.1.5 READOUT DEVICE.




FIG 6.Components of spectrophotometer.

2.1.1 SOURCE OF LIGHT

- Part of the UV and Visible radiation source is **Tungsten lamp**. See figure 7.



FIG.7.Tungsten lamp

- UV radiation source is **Deuterium or Hydrogen lamp** . See figure 8.
- Range of wavelength 200-400 nm.



FIG.8.Deuterium lamp

2.1.2 MONOCHROMATOR

- It is a device that breaks the polychromatic radiation into component wavelengths. See figure 9.




FIG.9.Monochromator components.

The monochromator unit consists of :

- **Entrance slit:** defines narrow beam of radiation from source.
- **Collimating mirror:**(polished surface) collimates the lights.
- **Diffraction grating or Prism** (make of quartz): disperses the light into specific wavelength.
- **Focusing mirror:** captures the dispersed light & sharpens the same to the sample via exit slit
- **Exit slit:** allows the corrected wavelength of light to the sample .



2.1.3 SMPLE SOLIOTION IN CUVETTE

- liquid sample is usually contained in a cell called a *cuvette*. See figure10.
- Fingerprints and droplets of water disrupt light rays during measurement.
- Cuvette from *Quartz* can be used in UV as well as in visible spectroscopy.
- Cuvette from *Glass* is suitable for visible but not for UV spectroscopy because it absorbs UV radiation.



FIG 10.sample solution in cuvette

2.1.4 PHOTO DETECTOR

- A photo detector is a semiconductor device which converts light energy to electrical energy. It consists of a simple P-N junction diode and is designed to work in reverse biased condition. The photons approaching the diode are absorbed by the photodiode and current is generated.⁴ See figure 11.




FIG 11. Photodiode

2.1.5 READOUT DEVICE.

- Digital screen to record an uv spectrograph with absorbance against the wavelength.

2.2 TYPES of SPECTROPHOTOMETER




FIG 12. Types of spectrophotometer.

2.3 APPLICATIONS OF UV-VIS SPECTROSCOPY.

- UV-Vis spectroscopy is used heavily in many different research areas to identify or quantify a sample.³
- Chemical field:
 1. Detection of impurities.
 2. Structure of organic compounds (single or double bond, presence or absence of functional group).
 3. Kinetics of reaction.
 4. Manufacturing drugs.
- Biological fields
 1. quantify the amount of protein and DNA in a sample
 2. quantify the amount of bacterial cells in a cell culture

○ Major **advantages** of uv-vis spectroscopy are:

1. High sensitivity.
2. Require only small volume of sample.
3. Linearity over wide range of concentration.
4. Can be used with gradient elution.⁴

○ Major **disadvantages** of uv-vis spectroscopy are:

1. Not linear for high concentration.
2. Does not work with compounds that do not absorb light at this wavelength region.
3. Generates significant heat and requires external cooling.⁴



○ Thank you

